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
5-B

Diels-Alder Reaction

And

Some Related Reactions

Diels-Alder Reaction

Reagent for reaction 1

Maleic anhydride (mw=98.06) (2.0 g, 20.2 mmol)
Cyclopentadiene (mw=66.10) (2.0 mL, 24.2 mmol)

AcOEt (7 mL)
Ligroin (7 mL)

Preparation of reaction 1

1. Maleic anhydride and ethyl acetate (AcOEt) were placed in a 50 mL Erlenmeyer flask. The mixture was swirled on a hot plate at 55 °C until the maleic anhydride was completely dissolved.
2. The solution was then cooled in an ice bath. After cooling, cyclopentadiene was added, and the mixture was swirled for over 15 minutes at that temperature.
3. After that, the mixture was heated with a hot plate at 90 °C until the resulting precipitate was completely dissolved.
4. Subsequently, the solution was allowed to cool to room temperature, during which crystals formed and were collected by filtration.

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
Product 1 (mw=164.16) (1.5 g, 9.1 mmol)
water (20 mL)

Preparation of reaction 2

- Product 1 and 20 mL of water were placed in a 100 mL Erlenmeyer flask. The mixture was stirred with a magnetic stir bar on a hot plate set to 120 °C for 30 min or longer.
- After that, the reaction mixture was allowed to cool to room temperature until crystals began to form. Once crystallization started, the mixture was further cooled in an ice bath to promote crystal growth.
- The resulting crystals were collected by filtration.


Regent for reaction 3

Product 2 (mw=182.18) (1.0g, 9.1 mmol)
conc.H₂SO₄ (5 mL)

Preparation of reaction 3

- Product 2 and conc. H₂SO₄ were placed in 50 mL Erlenmeyer flask. The mixture was swirled on a hot plate at 120 °C until all precipitate was dissolved.
- The mixture was cooled in an ice bath. After that, ice water was added to the mixture.
- The mixture was heated on a hot plate at 120 °C until all precipitate was dissolved.
- After that, the mixture was cooled in an ice bath.
- Resulting precipitate was filtrate. If crystallization does not occur or is insufficient, the inner wall of the Erlenmeyer flask may be gently scratched with a spatula to induce crystallization.

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